



Designed specifically for lumbar and cervical traction procedures, this model features a rolling top to maximise the effectiveness of applied traction forces for both static and intermittent traction procedures. The table is available with either electric or hydraulic height operation. A platform with adjustable mast is fitted for the addition of the optional traction machine and suitable for use with several makes of traction unit. A cervical in sitting unit is also supplied as an alternative to traditional cervical traction with the patient in a supine position.




Standard Features

- UK manufactured
- Lifetime frame & 5 year parts guarantee
- Choice of electric or hydraulic height operation
- Safe working load 250kg
- Adjustable angle backrest from horizontal to +85° with fittings for cervical in sitting procedures
- Positive and negative adjustment of head section from -20° to +85°
- Adjustable height traction machine mount included
- Cervical in sitting accessory, traction mast and all harnesses included
- Retractable wheel system
- Flexion stool optional with table but included with the Traction Package
- Low patient transfer height of 47cm (18")
- Wide choice of options, colours, and accessories

Overall Dimensions

Length: 198cm (78") approx
 Width: 63cm (25") approx
 Height: 47cm - 97cm approx

Upholstery Dimensions



Options & Accessories




Flexion Stool



Paper Roll Holder



Traction Machine



Battery Backup




Footswitch



Spare Harness Set

Operational Parameters



EMS Physio Ltd order codes
COU005 Hydraulic with breathing slot
COU006 Electric with breathing slot

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note; Unless specified Sapphire upholstery will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond

For Your Piece of Mind

- Lifetime frame and 5 year parts (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average life span of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2015 by LQRA

Standards

- IP Rating - IPX6
- EE Classification - EU MDR 2017/745 Annex VIII Class 1
- Electrical Standard - IEC 60601-1
- Vinyl flame retardant to CRIB 5 (BS 5852)
- Foam flame retardant to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required.

EMS PHYSIO Ltd
Downsview Road
Wantage
Oxfordshire OX12 9FE

t: +44 (0) 1235 772272
f: +44 (0) 1235 763518
e: sales@emsphysio.co.uk
www.emsphysio.co.uk